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S&T®
Basic
informations
*Grund-
Information*

The Perfect Hole

A good arteriotomy is essential for an end-to-side anastomosis. The hole should be elliptical and clean-cut, with even, steeply sloping edges. The hole should not have an irregular, ragged edge as it may weaken the wall and attract thrombosis. A well-made hole allows the surgeon to see the edges clearly, place the needle accurately, and tie sutures securely. Such a hole is hard to make, as every microvascular surgeon knows. Scissors just won't cut it.

Why Scissors Won't Cut it

The thick wall of a big vessel is tough and slippery. It has two layers which can separate with unappreciated ease, especially in older patients. Scissors don't just cut the vessel wall — they squeeze it, shear it, push it away, then cut it. Their action is unpredictable. Sometimes it's alarming. Sometimes it's dangerous. To make an accurate cut with scissors is difficult. Worst of all, scissors can separate a vessel wall into its two layers, cutting only the outer one and squeezing the inner one, uncut, into the lumen.

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Arteriotomy Clamp

Prof. Dr. Robert D. Acland, University of Louisville, Louisville, Ky, USA
Dr. Joseph Banis, University of Louisville, Louisville, Ky, USA

The Acland-Banis Arteriotomy Clamp was developed to overcome the difficulties of creating a hole in the side of a vessel, in an end-to-side anastomosis. The device enables the surgeon to make an arteriotomy with an instrument that cuts the vessel the best — a knife. The arteriotomy clamp picks up and holds the piece of the vessel wall that is to be excised. Then, with a single curving sweep of a straight blade that is held close against the clamp tip, the ellipse is excised. The result is a perfect arteriotomy, every time.

The clamp is a precision-made, cross-action forceps with a heavy closing pressure and a tapered, rounded, fenestrated tip. The hole in the tip is important because the tissue that is picked up can squeeze into it and be held firmly. There are three sizes of clamps, with tips 1.2, 1.8 and 2.4 mm in diameter. Each clamp can cut a different size hole, depending on the depth of the tissue being picked up. Using the set of three clamps, holes can be cut with mean diameter ranging from 1 to 5 mm.

Arteriotomieklemme

Prof. Dr. Robert D. Acland, University of Louisville, Louisville, Ky, USA
Dr. Joseph Banis, University of Louisville, Louisville, Ky, USA

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
Basic informations Grund- Information

Using the Clamp

1. Mobilize the vessel so that it can be rotated freely.
2. Apply occluding clamps, keeping the vessel at its natural state of resting tension.
3. Empty the occluded segment by finger pressure while momentarily releasing one clamp
4. Remove adventitia thoroughly and widely around the planned arteriotomy site.
5. Rotate the vessel so that the site of arteriotomy is seen in clear profile.
6. Have the small vessel in view, with its end prepared and dilated, in order to judge its size.
7. With an arteriotomy clamp of appropriate size, carefully pick up a piece of vessel wall which is half as long and just as wide as the desired arteriotomy.
Note that because of the pressure which the tip of the clamp exerts, this step is irreversible.
8. With a CBS-35 blade held in metal-to-metal contact with the clamp tip, cut around the tip with a steady sweeping action.

Important!

With a practiced use, the size of arteriotomy which results is highly predictable. However, we strongly recommend a brief initial practice period on animal or fresh cadaver vessels. These vessels should be stretched out to their natural length and tension to simulate real conditions. Remember that the arteriotomy will be twice as long and just as wide as the piece which is picked up and excised.

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Arteriotomy Clamp by Acland-Banis, Model registered

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

S&T[®]
Basic
informations
*Grund-
Information*

**Genuine TUPPERHAND
Universal Hand Holder
and Retractor Set**

Jack Tupper, M.D., Oakland, CA, USA

The hand surgeon's best friend:

Tupper's Universal Hand Holder and Retractor Set enables the surgeon to fixate the hand in any position needed. Ball chains of varying lengths, with skin or tendon hooks, fixate the hand securely in the handholder. With the patient's hand held firmly, the surgeon can work unhindered beneath a microscope.

The slots around the handholder facilitate retraction in nearly every direction. Upward retraction is easily achieved with the elevated chain holders. The special container is designed for easy handling and safe storage of the set in and outside of the O.R. All parts of the **TUPPERHAND** are made of stainless steel, cleaning and sterilisation procedures are the same as with other surgical instruments.

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


Best protection

The TH 110 container was specifically designed to store the Tupper handtable and accessories. Made of the highest quality stainless steel, these cases afford your instruments the best protection possible.

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00609

S&T[®]
Basic
informations
*Grund-
Information*

Surgery of Peripheral Nerves and Tendons

by Meyer

Prof. Dr. Viktor E. Meyer, University Hospital Zurich, Zurich, Switzerland

In reconstruction of peripheral nerves, the exact fascicular coaptation is one of several important factors, upon which we have a direct influence, at least to a certain degree.

A straight and even cut-surface of the nerve stumps greatly facilitates the exact fascicular coaptation.

The submitted holding forceps were developed to trim nerve stumps smoothly with a blade, without additional trauma by the holding device and with a minimal loss of tissue. The simple design makes their application possible also in cases of difficult surgical exposure. The eccentrically placed guide slot makes a very economical trimming of the nerves possible, for example, as preparation for a direct suture.

Scanning electron microscopical comparison of stumps of the human sural nerve, cut by a toothed pair of scissors, respectively by the method applied here, have shown that the anatomical changes of the cut-surface is less evident by using a sharp blade than by using a toothed pair of scissors (V. E. Meyer and J. Smahel: The Surgical Cut-Surface of Peripheral Nerves, International Journal for Microsurgery, 2: 187, 1980). By this technique smooth and straight surface-cuts can be achieved, especially on larger nerves, which also facilitates assessment of the fascicular topography.

For trimming of small single fascicles or small groups of fascicles the application of serrated scissors is recommended.

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Basic informations

Grund-Information


Figure 1:
Scanning electron microscopical picture of a resection by above technique of a non-fixed human flexor tendon.

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
Figure 2:
The holding forceps is a gentle handling instrument for the tendon stumps at suturing.

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
Figure 3:
Resected slices of median nerve at a secondary suturing. Using the holding forceps slices as thin as 1 mm can easily be cut.

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Clinical and laboratory experience has shown that this instrumentarium facilitates the dissection and/or the trimming of other cylindrical anatomical structures, above all, if in a saving manner a smooth surface should be achieved (Tendons, Fallopian Tubes, Vas Deferens, Vessels, etc.): **It has proven its usefulness especially in hand surgery at trimming the tendon stumps.** Not only extraordinarily smooth tendon surfaces can be cut, but the instrumentarium also permits the gentle handling of the stumps for suturing (Fig. 1 and Fig. 2).

Applied holding forceps, with blade in guide slot. An oscillating movement under slight pressure and constant rinsing results in a straight, smooth resection of the nerve end.

Insert blade carefully into guide slot to avoid damage to its cutting edge.

Clean the guide slot under running water, with the back of the used blade. Avoid spreading and distorting the guide slot. Don't use blades, which are too thick.

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

The examination of cut-surfaces of peripheral nerves and tendons has been performed in cooperation with J. Smahel, M.D., Chief of the Research Laboratory of the Surgical Clinic B, University of Zurich Medical School, Zurich, Switzerland.

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

Nerve and Tendon Sets, Instrument Rack

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
"The special Vascular Dissecting Forceps of **S&T** are used as fine dissectors for precision work around the small vessels such as perforators during operation involving perforator flaps and other flaps. These dissectors allow for precision dissection with minimal trauma to the adjacent tissue."

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Delicate tips and blades
Feine Spitzen und Klingen

Accessories Zubehör

Acland's Practice Manual For Microvascular Surgery

By Robert D. Acland and S. Raja Sabapathy

The benchmark for everybody who is looking for a high level of basic skills in microvascular surgery.

This practice manual includes amongst others: basic instrumentation, microscope adjustment, suturing, rat anaesthesia, arterial and venous anastomosis, interpositional vein graft, one-way-up and continuous suture anastomosis.

New and revised edition of 2008.

Important!
Exclusively available at your **S&T**[®] distributor.

Acland's Practice Manual For Microvascular Surgery

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
Acland's Practice Manual For Microvascular Surgery, 3rd edition
Acland's Practice Manual For Microvascular Surgery, 3. Auflage

99953