


Övriga instrument | Various instruments | Sonstige
Instrumente | Altri strumenti | Instrumentos varios |
様々な器具

Ergonomic design ensures
comfortable grip and exact
handling of delicate tissue.

All instruments are manufactured in
selected top quality surgical steel.





Gillquist-STILLE artroskopisugrör | Gillquist-STILLE Arthroscopy Suction Tube

Arthroscopy suction tube with mandrin.

119-10915 28 cm 11"

Oretorp-STILLE sond | Oretorp-STILLE Probe

Blunt ended probe for arthroscopy. Graduated shaft.

119-10928 Blunt 19.5 cm 7.75"



STILLE sond | STILLE Probe

Blunt ended probe for arthroscopy. Graduated shaft.

119-10928-1 Blunt 19.5 cm 7.75"



STILLE uterusond | STILLE Uterine Sound

Uterine sound with malleable silver shaft. Graduated in cm.
Drop-shaped tip.

119-28452

27 cm 10.75"



Bergström-STILLE muskelbiopsi kanyl | Bergström-STILLE Muscle biopsy Cannula

Muscle biopsy cannula. Consists of a sharptipped, hollow outer needle with a small opening near the tip. A punch with a sharp edge fits tightly into the needle. The stylet is used for pushing out biopsy specimens from the punch.

119-29187-50 Diam 5 mm

119-29187-60 Diam 6 mm



STILLE Jameson skjutmått | STILLE Jameson Caliper

Caliper graduated 0-80 mm.

119-30025 0-80 mm 0-3" 9.5 cm 3.75"



STILLE Castroviejo skjutmått | STILLE Castroviejo Caliper

Caliper with locking screw, graduated 0-20 mm.

119-30022 0-20 mm 8 cm 3.1"



STILLE mätsticka | STILLE Ruler

Stainless steel ruler graduated in mm, 0-15 cm.

119-30010-15 15 cm 6"





STILLE skalpellhandtag | STILLE Scalpel Handle

The most common scalpel handles made with highest precision.

119-50010-3	No 3	Blade 10-15	12 cm 5"
119-50010-13	No 3L	Blade 10-15	21 cm 8.25"
119-50010-4	No 4	Blade 18-25	13 cm 5.25"
119-50010-14	No 4L	Blade 18-25	21 cm 8.25"
119-50010-7	No 7	Blade 10-15	16 cm 6.25"



STILLE Barron skalpellhandtag | STILLE Barron Scalpel Handle

8-sided handle for precise handling and rotation.

119-50012	Blade 10-15	14 cm 5.5"
-----------	-------------	------------

STILLE skalpellhandtag | STILLE Scalpel Handle

Round handle guarantees precise finger tip control.


119-50015	Blade 10-15	14.5 cm 5.75"
-----------	-------------	---------------



STILLE Mayo peanghängare | STILLE Mayo Safety Pin

Safety Pin with ball for better control of instruments on the tray.

119-50022 13.5 cm 5.25"



STILLE Hemostop | STILLE Hemostop

Designed as an aid to stemming complicated venous and arterial bleeding during open surgery. Relatively long to be able to reach down into the pelvis.

119-60001 Ring diameter 18/28 mm 24 cm 9.5"





Medicway instrumentolja | Medicway Instrument Oil

Instrument oil, package of 2 pcs à 50 ml.

119-78900-1 50 ml, package of 2 pcs

STILLE rengöringsborste | STILLE Cleaning Brush

Cleaning brush for surgical instruments. Special stainless steel fine threaded, bristle diameter 0.12 mm, height 20 mm. Handy grip of plastic material. Autoclavable at 121C.


119-78950 20 cm 8", package 2 pcs



STILLE Instrumentbricka | STILLE Instrument Sterilization Tray

Instrument tray for transport, sterilization and storage of delicate instruments. The trays are designed to assure proper penetration of sterilizing agent between the mat and the tray surface. Silicone mat included. Available in three different sizes. The trays are made to withstand common sterilization methods, including steam, gas, cold sterilization solutions, and temperatures up to 182° C/360 F.


199-60010	Micro tray	10.1x19x1.9 cm	4x7.5x0.75"
199-60020	Micro tray	15.2x25.4x1.9 cm	6x10x0.75"
199-60030	Instrument tray	25.4x38.1x1.9 cm	10x15x0.75"



STILLE-Gigli trådsåg | STILLE-Gigli Wire Saw

Gigli wire saw with superior cutting quality. Smooth joints between loops and wire. Available in three different lengths.

125-19150-30	Pack of 25 pcs	30 cm 12"
125-19150-40	Pack of 25 pcs	40 cm 15.75"
125-19150-50	Pack of 25 pcs	50 cm 20"
125-19160	Handle	6.5 cm 2.5"
125-19170	Guide/Dura Protector	34 cm 13"



STILLE gipssåg | STILLE Plaster Cast Remover

A complete kit of plaster cast removers including wire saws in two lengths, stainless steel handles and plastic profiles.

413-2500-2	Metal Handle		
413-2502-85	Wire Saws	100 pcs	85 cm 33.5"
413-2502-125	Wire Saws	100 pcs	125 cm 49.25"
413-2504-85	Soft Profiles	100 pcs	85 cm 33.5"
413-2504-125	Soft Profiles	100 pcs	125 cm 49.25"



STILLE frakturtråd | STILLE Bone Wire

STILLE bone wire on 5 m coils. Not available in USA.

121-18013-4	Stainless	0.4 mm	5 m coils diam ,soft
121-18013-5	Stainless	0.5 mm	5 m coils diam, soft
121-18013-6	Stainless	0.6 mm	5 m coils diam, soft
121-18013-7	Stainless	0.7 mm	5 m coils diam, soft
121-18013-8	Stainless	0.8 mm	5 m coils diam, soft
121-18013-9	Stainless	0.9 mm	5 m coils diam, soft
121-18013-10	Stainless	1.0 mm	5 m coils diam, soft



Palmer-STILLE benstift | Palmer-STILLE Bone Nail

Round bone nails with diameter 1.5 mm. Not available in USA.

121-18224-15	Stainless	Round	Length 15 mm	Pkg of 5 pcs
121-18224-20	Stainless	Round	Length 20 mm	Pkg of 5 pcs
121-18224-25	Stainless	Round	Length 25 mm	Pkg of 5 pcs
121-18224-30	Stainless	Round	Length 30 mm	Pkg of 5 pcs
121-18224-35	Stainless	Round	Length 35 mm	Pkg of 5 pcs
121-18224-40	Stainless	Round	Length 40 mm	Pkg of 5 pcs
121-18224-45	Stainless	Round	Length 45 mm	Pkg of 5 pcs
121-18224-50	Stainless	Round	Length 50 mm	Pkg of 5 pcs
121-18224-55	Stainless	Round	Length 55 mm	Pkg of 5 pcs
121-18224-60	Stainless	Round	Length 60 mm	Pkg of 5 pcs
121-18224-65	Stainless	Round	Length 65 mm	Pkg of 5 pcs
121-18224-70	Stainless	Round	Length 70 mm	Pkg of 5 pcs
121-18227	Holder/Driver		12 cm 5.7"	
121-18229	Impactor		17 cm 6.75"	





Wiberg-STILLE benmärkla | Wiberg-STILLE Fracture Staple

Fracture staples in three different sizes. Not available in USA.

121-18251-16	Stainless	16 mm	Pkg of 5 pcs
121-18251-19	Stainless	19 mm	Pkg of 5 pcs
121-18251-22	Stainless	22 mm	Pkg of 5 pcs
121-18255	Holder/Driver	17 cm 6.75"	

STILLE Kirschner trådborr | STILLE Kirschner Wire Drill

Wire drills with lancet points with or without eye. Not available in USA.

121-19357-11	Without eye	1.1 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-13	Without eye	1.3 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-15	Without eye	1.5 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-18	Without eye	1.8 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-20	Without eye	2.0 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-213	With eye	1.3 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-215	With eye	1.5 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-218	With eye	1.8 mm	23.5 cm 9.25"	Pkg of 5 pcs
121-19357-220	With eye	2.0 mm	23.5 cm 9.25"	Pkg of 5 pcs

