


Slevar, Raspatorier, Elevatorier etc. | Currettes,  
Raspatories, Elevators etc. | Löffel, Raspatorien,  
Elevatorien etc. | Curette, Staccaperiostio,  
Elevatori etc. | Cucharillas, Legras, Elevadores etc. |  
鋭匙、やすり 骨膜起子




Tips stay sharp  
for a long time.

High HRC ensures strong  
and sharp blades.

Handles designed to  
prevent from slipping.

## Récamier STILLE uteruskyrett | Récamier STILLE Uterine curette

Malleable uterine curettes in different sizes. With shaft in non-hardened stainless steel, ensures long-time function without any surface coating falling off. Higher patient safety and longer life-time expectancy.



111-20340-5

111-20340-7

111-20340-8


111-20340-10



111-20350-7

111-20350-8

111-20350-10



111-20350-12


111-20340-5	Sharp	5 mm	30 cm 12"
111-20340-7	Sharp	7 mm	30 cm 12"
111-20340-8	Sharp	8 mm	30 cm 12"
111-20340-10	Sharp	10 mm	30 cm 12"
111-20350-7	Blunt	7 mm	30 cm 12"
111-20350-8	Blunt	8 mm	30 cm 12"
111-20350-10	Blunt	10 mm	30 cm 12"
111-20350-12	Blunt	12 mm	30 cm 12"





### STILLE slev | STILLE Bone curette

Oval and sharp bone curettes with widths from 2 mm to 12 mm.




111-20502-2	Sharp oval	2 mm	17,5 cm 7"
111-20502-3	Sharp oval	3 mm	18 cm 7"
111-20502-4	Sharp oval	4 mm	18,5 cm 7.25"
111-20502-5	Sharp oval	5 mm	18,5 cm 7.25"
111-20502-6	Sharp oval	6 mm	18,5 cm 7.25"
111-20502-7	Sharp oval	7 mm	19 cm 7.5"
111-20502-8	Sharp oval	8 mm	19,5 cm 7.75"
111-20502-10	Sharp oval	10 mm	21 cm 8.25"
111-20502-12	Sharp oval	12 mm	21 cm 8.25"



### STILLE slev | STILLE Bone curette

Oval and sharp bone curettes with widths from 2 mm to 6 mm. Longer Shaft version.



111-20503-2	Sharp oval	2 mm	23 cm 9"
111-20503-3	Sharp oval	3 mm	23 cm 9"
111-20503-4	Sharp oval	4 mm	23 cm 9"
111-20503-5	Sharp oval	5 mm	23 cm 9"
111-20503-6	Sharp oval	6 mm	23 cm 9"



### Olivecrona STILLE raspatorium | Olivecrona STILLE Raspatory

Raspatory with sturdy and ergonomic handle. Curved tip.



111-22150

111-22150	Curved	11 mm	17 cm 6.75"
-----------	--------	-------	-------------

### Cobb STILLE raspatorium | Cobb STILLE Raspatory

Popular raspatory with round handle. Available with several different tip sizes.

111-22155-10	10 mm	28,5 cm 11.25"
111-22155-13	13 mm	28,5 cm 11.25"
111-22155-19	19 mm	28,5 cm 11.25"
111-22155-25	25 mm	28,5 cm 11.25"





### Cottle STILLE raspatorium | Cottle STILLE Raspatory

Commonly used raspatory with a sharp 9 mm rounded tip.

111-22158-9 9 mm 19 cm 7.5"



### Lambotte STILLE raspatorium | Lambotte STILLE Raspatory

Raspatory with a wide 23 mm tip.




111-22160

111-22160 Curved 23 mm 21 cm 8.25"



### Wiberg STILLE raspatorium | Wiberg STILLE Raspatory

A straight raspatorium in four different sizes.



111-22171-4

111-22171-7

111-22171-9

111-22171-12

111-22171-4	Straight	4 mm	17 cm 6.75"
111-22171-7	Straight	7 mm	17 cm 6.75"
111-22171-9	Straight	9 mm	17 cm 6.75"
111-22171-12	Straight	12mm	17 cm 6.75"



### Ragnell STILLE raspatorium | Ragnell STILLE Raspatory


Double-ended raspatorium with blunt hooks.

111-22350 Blunt hook 17 cm 6.75"



### Howarth STILLE näsraspatorium | Howarth STILLE Nasal Raspatory

Double-ended nasal raspatory with sharp hooks.



111-22370

111-22370 Double-ended 22 cm 8.75"




### Joseph STILLE raspatorium | Joseph STILLE Raspatory

Slightly curved raspatory with flat handle.

111-22385-4 4 mm 17 cm 6.75"

### STILLE näsraspatorium | STILLE Nasal Raspatory

Sharp, slightly curved nasal raspatory.



111-22520-4

111-22520-4 Slightly curved Sharp 4 mm 18 cm 7"







**Freer STILLE elevatorium | Freer STILLE Elevator**


Double-ended elevator sharp/blunt with rounded, delicate handle.

111-22140-4	4.5 mm	19.5 cm 7.75"
111-22140-6	6.5 mm	19.5 cm 7.75"



**STILLE näselevatorium | STILLE Nasal Elevator**

Blunt nasal elevator.



111-22522-4

111-22522-4	Slightly curved blunt	4 mm	18 cm 7"
-------------	-----------------------	------	----------



**Cottle STILLE septum elevator | Cottle STILLE Septum Elevator**

Double-ended fine elevator with graduated shafts.

111-22532	Graduated	3.3/3.7 mm	23 cm 9"
-----------	-----------	------------	----------



**Langenbeck STILLE elevator | Langenbeck STILLE Elevator**

Sturdy elevator mainly used in orthopaedics. Available in two sizes.

111-22550-8	8 mm	19.5 cm 7.75"
111-22550-10	10 mm	19.5 cm 7.75"



### Mini-Hohmann STILLE elevatorium | Mini-Hohmann STILLE Elevator

Finer version of the Hohmann elevator. Used in foot surgery.

111-22565-6	Bone lever	6 mm	16 cm 6.25"
111-22565-8	Bone lever	8 mm	16 cm 6.25"



### Hohmann STILLE elevatorium | Hohmann STILLE Elevator

General utility bone lever for orthopaedic surgery.

111-22567-9	Rounded	9 mm	22 cm 8.75"
111-22568-18	Rounded	18 mm	24 cm 9.5"
111-22569-18	Flat	18 mm	25,5 cm 10"
111-22569-24	Flat	24 mm	26 cm 10.25"



### STILLE elevatorium | STILLE Elevator

Standard bone lever in two different sizes.

111-22563-3	3 mm	14,5 cm 5.75"
111-22563-6	6 mm	14,5 cm 5.75"



### Stolte STILLE tonsillelevatorium | Stolte STILLE Tonsil Elevator

Curved double-ended tonsil elevator.




111-22635

111-22635	Double-ended	20 cm 8"
-----------	--------------	----------

## STILLE dissektor | STILLE Dissector

Double-ended dissector in several different models.



111-22440-18	Fine with hole	18 cm 7"
111-22440-24	Fine with hole	24 cm 9.5"
111-22442-19	Coarse with hole	19.5 cm 7.75"
111-22442-24	Coarse with hole	24 cm 9.5"
111-22444	Fine without hole	24.5 cm 9.75"
111-22446	Coarse without hole	24.5 cm 9.5"



## Fomon STILLE rasp | Fomon STILLE Rasp

Double-ended rasp with TC-inserts for better durability. Low profile for better access in narrow areas.

111-22860	TC-inserts 1&2	21.5 cm 8.5"
111-22861	TC-inserts 3&4	21.5 cm 8.5"
111-22862	TC-inserts 5&6	21.5 cm 8.5"
111-22863	TC-inserts 7&8	21.5 cm 8.5"

