


Nålförare | Needle holders | Nadelhalter | Portaghi | Porta-agujas | 持針器





Nålförare storlekstabell | Needle holder in relation to suture size | Nadelhalter in Relation zur Stärke des Nahtmaterials |
 Il portaaghi in relazione alla misura della sutura | Porta-agujas en relación al tamaño de la sutura | 縫合糸サイズに適する持針器



Art.no	Description	Max suture size
104-15542	RAPP-STILLE KNOT-HOLDING FORCEPS 13 cm 5"	3/0
104-15545	WEBSTER STILLE very delicate DD jaws 12 cm 4.75"	5/0
104-15546	WEBSTER STILLE very delicate smooth 12 cm 4.75"	5/0
104-15549-131	KILNER STILLE TC jaws 13 cm 5"	4/0
104-15550	GILLIES STILLE serrated 16 cm 6.25"	2/0
104-15551	GILLIES STILLE TC jaws 16 cm 6.25"	2/0
104-15556	WEBSTER STILLE delicate smooth 12 cm 4.75"	4/0
104-15557	WEBSTER-KLEINERT STILLE smooth offset bows 12 cm 4.75"	4/0
104-15558	WEBSTER STILLE smooth TC jaws 12 cm 4.75"	3/0
104-15559	WEBSTER STILLE serrated TC jaws 12 cm 4.75"	4/0
104-15560	HALSEY STILLE TC jaws 12.5 cm 5"	3/0
104-15560-2	HALSEY STILLE ERGONOMIC TC jaws 12.5 cm 5"	3/0
104-15561-141	MAYO-HEGAR STILLE TC jaws 14 cm 5.5"	0
104-15561-161	MAYO-HEGAR STILLE TC jaws 16 cm 6.25"	0
104-15561-181	MAYO-HEGAR STILLE TC jaws 18 cm 7"	0
104-15561-201	MAYO-HEGAR STILLE TC jaws 20 cm 8"	0
104-15561-251	MAYO-HEGAR STILLE TC jaws 25 cm 10"	0
104-15562-161	MAYO-HEGAR STILLE Heavy Pattern TC jaws 16 cm 6.25"	1
104-15562-181	MAYO-HEGAR STILLE Heavy Pattern TC jaws 18 cm 7"	1
104-15562-201	MAYO-HEGAR STILLE Heavy Pattern TC jaws 20 cm 8"	1
104-15562-241	MAYO-HEGAR STILLE Heavy Pattern TC jaws 24 cm 9.5"	1
104-15562-271	MAYO-HEGAR STILLE Heavy Pattern TC jaws 27 cm 10.5"	1
104-15562-301	MAYO-HEGAR STILLE Heavy Pattern TC jaws 30 cm 12"	1
104-15565-151	CRILE-WOOD STILLE TC jaws 15 cm 6"	2/0
104-15565-181	CRILE-WOOD STILLE TC jaws 18 cm 7"	2/0
104-15565-201	CRILE-WOOD STILLE TC jaws 20 cm 8"	2/0
104-15565-231	CRILE-WOOD STILLE TC jaws 23 cm 9"	2/0

Art.no	Description	Max suture size
104-15566-151	BABY CRILE-WOOD STILLE delicate TC 15 cm 6"	3/0
104-15566-181	CRILE-WOOD STILLE delicate TC jaws 18 cm 7"	3/0
104-15566-201	CRILE-WOOD STILLE delicate TC jaws 20 cm 8"	3/0
104-15566-231	CRILE-WOOD STILLE delicate TC jaws 23 cm 9"	3/0
104-15571	NEW ORLEANS STILLE TC jaws 18 cm 7"	3/0
104-15573-161	RYDER STILLE TC jaws 16 cm 6.25"	3/0
104-15573-181	RYDER STILLE TC jaws 18 cm 7"	3/0
104-15573-201	RYDER STILLE TC jaws 20 cm 8"	3/0
104-15585-161	STILLE DD jaws 16 cm 6.25"	4/0
104-15585-181	STILLE DD jaws 18.5 cm 7.5"	4/0
104-15585-221	STILLE DD jaws 22 cm 8.75"	4/0
104-15585-251	STILLE DD jaws 25 cm 10"	4/0
104-15586	STILLE TC jaws 18 cm 7"	4/0
104-15590-121	OLSEN-HEGAR STILLE TC jaws with scissors 12 cm 4.75"	3/0
104-15590-131	OLSEN-HEGAR STILLE TC jaws with scissors 13.5 cm 5.25"	2/0
104-15590-171	OLSEN-HEGAR STILLE TC jaws with scissors 17 cm 6.75"	0
104-15590-181	OLSEN-HEGAR STILLE TC jaws with scissors 18.5 cm 7.25"	0
104-15621-179	MATHIEU STILLE TC jaws 17 cm 6.75"	1
104-15621-209	MATHIEU STILLE TC jaws 20 cm 8"	1
104-15665-151	Frimand-STILLE DUO-GRIP Needle holder TC jaws 15 cm 6"	2/0
104-15777	KAYE STILLE with scissors serrated TC jaws 12 cm 4.75"	5/0
104-15785	CONVERSE STILLE TC jaws 10.5 cm 4.25"	4/0

Suture size indicates the largest recommended for each needle holder.

Derf STILLE nålförare | Derf STILLE Needle holder

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.



104-15542

104-15542 TC jaws 12 cm 4.75"



Webster STILLE nålförare | Webster STILLE Needle holder

Very delicate needle holder for very fine sutures. Smooth jaws or DD jaws. Intended for 5/0 suture and finer.




104-15545



104-15546

104-15545 Very delicate DD jaws 12 cm 4.75"

104-15546 Very delicate Smooth 12 cm 4.75"



Kilner STILLE nålförare | Kilner STILLE Needle holder


Delicate needleholder with bayonet shaped shanks. TC jaws. Intended for 4/0 suture and finer.

104-15549-131 13 cm 5"




Gillies STILLE nålförare | Gillies STILLE Needle Holder

Needle holder with scissors. Ergonomic design on shanks and finger rings. One rather large ring handle is set at an angle for fitting the top of the thumb. Without ratchets. With or without TC jaws. The scissor part enables suturing and cutting with one instrument. Intended for 2/0 suture and finer.



104-15550




104-15551

104-15550	With scissors		16 cm 6.25"
104-15551	With scissors	TC jaws	16 cm 6.25"



Webster/Kleinert STILLE nålförare | Webster/Kleinert STILLE Needle holder

Delicate needle holder for fine needles. With smooth jaws or TC jaws. The Webster-Kleinert (104-15557) has offset bows for a comfortable grip. Intended for 4/0 suture and finer.




104-15556

104-15557



104-15559

104-15556	Delicate	Smooth	12 cm 4.75"
104-15557	Offset bows	Smooth	12 cm 4.75"
104-15559	Serrated	TC jaws	12 cm 4.75"



Webster STILLE nålförare | Webster STILLE Needle holder


Delicate needle holder for fine needles. With smooth TC jaws. Intended for 3/0 suture and finer.



104-15558

104-15558	Smooth	TC jaws	12 cm 4.75"
-----------	--------	---------	-------------





Halsey STILLE nålförare | Halsey STILLE Needle Holder

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb. For slightly larger needles than the Webster needle holder. Intended for 3/0 suture and finer.



104-15560

104-15560	TC jaws	12.5 cm 5"
104-15560-2 Ergonomic	TC jaws	12.5 cm 5"

Mayo-Hegar STILLE nålförare | Mayo-Hegar STILLE Needle holder

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.



104-15561


104-15561-141	TC jaws	14 cm 5.5"
104-15561-161	TC jaws	16 cm 6.25"
104-15561-181	TC jaws	18 cm 7"
104-15561-201	TC jaws	20 cm 8"
104-15561-251	TC jaws	25 cm 10"



Mayo-Hegar STILLE kraftig nålförare | Mayo-Hegar STILLE Heavy Needle holder

Heavier version of the widely used Mayo-Hegar needle holder, with TC jaws. Intended for 1 suture and finer.

104-15562-161	Heavy Pattern	TC jaws	16 cm 6.25"
104-15562-181	Heavy Pattern	TC jaws	18.5 cm 7.25"
104-15562-201	Heavy Pattern	TC jaws	20.5 cm 8"
104-15562-241	Heavy Pattern	TC jaws	24 cm 9.5"
104-15562-271	Heavy Pattern	TC jaws	27 cm 10.75"
104-15562-301	Heavy Pattern	TC jaws	30.5 cm 12"




STILLE Crile-Wood nålförare | STILLE Crile-Wood Needle holder

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.



104-15565

104-15565-151	TC jaws 15 cm 6 in
104-15565-181	TC jaws 18 cm 7 in
104-15565-201	TC jaws 20 cm 8 in
104-15565-231	TC jaws 23 cm 9 in



Crile-Wood STILLE gracil nålförare | Crile-Wood STILLE Delicate Needle holder

Delicate version of the popular Crile-Wood needle holder. With TC jaws. Intended for 3/0 suture and finer.



104-15566

104-15566-151	Delicate	TC jaws	15 cm 6"
104-15566-181	Delicate	TC jaws	18 cm 7"
104-15566-201	Delicate	TC jaws	20 cm 8"
104-15566-231	Delicate	TC jaws	23 cm 9"



STILLE-Frimand DUO-GRIP nålförare | STILLE-Frimand DUO-GRIP Needle holder

Designed in cooperation with Carl-Fredrik Frimand, M.D. Intended for precise use in the palmar grip without the need of grip change during suturing. This new and innovative needle holder reduces the suturing time and tissue manipulation as well as offering better ergonomics for the surgeon. Intended for 2/0 suture and finer.



104-15665-151

104-15665-151	TC jaws	15 cm 6"
---------------	---------	----------



New Orleans STILLE nålförare | New Orleans STILLE Needle holder

Short tapered TC jaws, one jaw straight serrated and the other cross-serrated.
Intended for 3/0 suture and finer.



104-15571

104-15571	TC jaws	18 cm 7"
-----------	---------	----------



Ryder STILLE nålförare | Ryder STILLE Needle holder

Needle holder with very narrow, non-tapered, serrated TC jaws. Popular in vascular surgery. Intended for 3/0 suture and finer.



104-15573

104-15573-161	TC jaws	16 cm 6.25"
104-15573-181	TC jaws	18 cm 7"
104-15573-201	TC jaws	20 cm 8"



STILLE nålförare | STILLE Needle holder

TC-sprayed delicate jaws. Intended for 4/0 suture and finer.



104-15585

104-15585-161	DD jaws	16 cm 6.25"
104-15585-181	DD jaws	18.5 cm 7.5"
104-15585-221	DD jaws	22 cm 8.75"
104-15585-251	DD jaws	25 cm 10"



STILLE nålförare | STILLE Needle holder

Delicate needle holder with TC jaws. Intended for 4/0 suture and finer.



104-15586

104-15586	TC jaws	18 cm 7"
-----------	---------	----------



Olsen-Hegar STILLE nålförare | Olsen-Hegar STILLE Needle holder

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer. See table "Needle holder in relation to suture size" for suitable suture size for each article code. 3+3 ratchets, 104-15590-171 has 2+2 ratchets.

104-15590-121	Scissors	12 cm 4.75"
104-15590-131	Scissors	13.5 cm 5.25"
104-15590-171	Scissors	17 cm 6.75"
104-15590-181	Scissors	18.5 cm 7.25"




Mathieu STILLE nålförare | Mathieu STILLE Needle holder

Needle holder with straight TC jaws. Plier type handle with ratchets that allows the surgeon to open and close the needle holder rapidly. Intended for 1 suture and finer.



104-15621

104-15621-179	TC jaws	17 cm 6.75"
104-15621-209	TC jaws	20 cm 8"



Kaye STILLE nålförare | Kaye STILLE Needle holder

Delicate needle holder with short, fine TC jaws and sharp scissors. The scissor part enables suturing and cutting with one instrument. Intended for 5/0 suture and finer.



104-15777

104-15777 Serrated TC jaws 12 cm 4.75"



Converse STILLE nålförare | Converse STILLE Needle holder

Short and delicate needle holder for fine work. With TC jaws. Intended for 4/0 suture and finer.

104-15785 TC jaws 10.5 cm 4.25"



Rapp-STILLE knuthållartång | Rapp-STILLE Knot-holding Forceps

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.



104-15140


104-15140 13 cm 5"



STILLE ståltrådsroterare | STILLE Wire Twisting Forceps

A good tool for twisting steel wires. With TC jaws.

104-15832	TC-inserts	15 cm 6"
104-15835	TC-inserts	18.5 cm 7.25"



STILLE ståltrådshållare | STILLE Wire Holding Forceps

For holding steel wires. Equipped with double-action joint and TC jaws. The jaws close parallel to secure a firm grip of the wire.

104-15854	Fine	18 cm 7"
104-15855	Medium	18 cm 7"
104-15856	Large	18 cm 7"





Nålförare storlekstabell | Needle holder in relation to suture size | Nadelhalter in Relation zur Stärke des Nahtmaterials |
 Il portaaghi in relazione alla misura della sutura | Porta-agujas en relación al tamaño de la sutura | 縫合糸サイズに適する持針器



Art.no	Description	Max suture size
104-15542	RAPP-STILLE KNOT-HOLDING FORCEPS 13 cm 5"	3/0
104-15545	WEBSTER STILLE very delicate DD jaws 12 cm 4.75"	5/0
104-15546	WEBSTER STILLE very delicate smooth 12 cm 4.75"	5/0
104-15549-131	KILNER STILLE TC jaws 13 cm 5"	4/0
104-15550	GILLIES STILLE serrated 16 cm 6.25"	2/0
104-15551	GILLIES STILLE TC jaws 16 cm 6.25"	2/0
104-15556	WEBSTER STILLE delicate smooth 12 cm 4.75"	4/0
104-15557	WEBSTER-KLEINERT STILLE smooth offset bows 12 cm 4.75"	4/0
104-15558	WEBSTER STILLE smooth TC jaws 12 cm 4.75"	3/0
104-15559	WEBSTER STILLE serrated TC jaws 12 cm 4.75"	4/0
104-15560	HALSEY STILLE TC jaws 12.5 cm 5"	3/0
104-15560-2	HALSEY STILLE ERGONOMIC TC jaws 12.5 cm 5"	3/0
104-15561-141	MAYO-HEGAR STILLE TC jaws 14 cm 5.5"	0
104-15561-161	MAYO-HEGAR STILLE TC jaws 16 cm 6.25"	0
104-15561-181	MAYO-HEGAR STILLE TC jaws 18 cm 7"	0
104-15561-201	MAYO-HEGAR STILLE TC jaws 20 cm 8"	0
104-15561-251	MAYO-HEGAR STILLE TC jaws 25 cm 10"	0
104-15562-161	MAYO-HEGAR STILLE Heavy Pattern TC jaws 16 cm 6.25"	1
104-15562-181	MAYO-HEGAR STILLE Heavy Pattern TC jaws 18 cm 7"	1
104-15562-201	MAYO-HEGAR STILLE Heavy Pattern TC jaws 20 cm 8"	1
104-15562-241	MAYO-HEGAR STILLE Heavy Pattern TC jaws 24 cm 9.5"	1
104-15562-271	MAYO-HEGAR STILLE Heavy Pattern TC jaws 27 cm 10.5"	1
104-15562-301	MAYO-HEGAR STILLE Heavy Pattern TC jaws 30 cm 12"	1
104-15565-151	CRILE-WOOD STILLE TC jaws 15 cm 6"	2/0
104-15565-181	CRILE-WOOD STILLE TC jaws 18 cm 7"	2/0
104-15565-201	CRILE-WOOD STILLE TC jaws 20 cm 8"	2/0
104-15565-231	CRILE-WOOD STILLE TC jaws 23 cm 9"	2/0

Art.no	Description	Max suture size
104-15566-151	BABY CRILE-WOOD STILLE delicate TC 15 cm 6"	3/0
104-15566-181	CRILE-WOOD STILLE delicate TC jaws 18 cm 7"	3/0
104-15566-201	CRILE-WOOD STILLE delicate TC jaws 20 cm 8"	3/0
104-15566-231	CRILE-WOOD STILLE delicate TC jaws 23 cm 9"	3/0
104-15571	NEW ORLEANS STILLE TC jaws 18 cm 7"	3/0
104-15573-161	RYDER STILLE TC jaws 16 cm 6.25"	3/0
104-15573-181	RYDER STILLE TC jaws 18 cm 7"	3/0
104-15573-201	RYDER STILLE TC jaws 20 cm 8"	3/0
104-15585-161	STILLE DD jaws 16 cm 6.25"	4/0
104-15585-181	STILLE DD jaws 18.5 cm 7.5"	4/0
104-15585-221	STILLE DD jaws 22 cm 8.75"	4/0
104-15585-251	STILLE DD jaws 25 cm 10"	4/0
104-15586	STILLE TC jaws 18 cm 7"	4/0
104-15590-121	OLSEN-HEGAR STILLE TC jaws with scissors 12 cm 4.75"	0
104-15590-131	OLSEN-HEGAR STILLE TC jaws with scissors 13.5 cm 5.25"	0
104-15590-171	OLSEN-HEGAR STILLE TC jaws with scissors 17 cm 6.75"	0
104-15590-181	OLSEN-HEGAR STILLE TC jaws with scissors 18.5 cm 7.25"	0
104-15621-179	MATHIEU STILLE TC jaws 17 cm 6.75"	1
104-15621-209	MATHIEU STILLE TC jaws 18.5 cm 7.25"	1
104-15777	KAYE STILLE with scissors serrated TC jaws 12 cm 4.75"	5/0
104-15785	CONVERSE STILLE TC jaws 10.5 cm 4.25"	4/0

Suture size indicates the largest recommended for each needle holder.