

Shaped finger rings for balance and comfort

Slim and flexible shanks for easy handling


Effortless and safe locking and unlocking with ratchets

Beveled edges to prevent suture from getting caught

Smooth and effortless action by calibrated box locks. No sharp edges.

Uniform, rounded tips

Durable jaws that are fitted with TC-inserts



ART. NUMBER	DESCRIPTION	MAX SUTURE SIZE	ART. NUMBER	DESCRIPTION	MAX SUTURE SIZE
104-15542	Derf STILLE Knot-holding Forceps 12 cm 4.75"	3/0	104-15566-151	Baby Crile-Wood STILLE Delicate TC 15 cm 6"	3/0
104-15545	Webster STILLE Very delicate Diamond jaws 12 cm 4.75"	5/0	104-15566-181	Crile-Wood STILLE Delicate TC jaws 18 cm 7"	3/0
104-15546	Webster STILLE Very delicate Smooth 12 cm 4.75"	5/0	104-15566-201	Crile-Wood STILLE Delicate TC jaws 20 cm 8"	3/0
104-15549-131	Kilner STILLE TC jaws 13 cm 5"	4/0	104-15566-231	Crile-Wood STILLE Delicate TC jaws 23 cm 9"	3/0
104-15550	Gillies STILLE Serrated 16 cm 6.25"	2/0	104-15571	New Orleans STILLE TC jaws 18 cm 7"	3/0
104-15551	Gillies STILLE TC jaws 16 cm 6.25"	2/0	104-15573-161	Ryder STILLE TC jaws 16 cm 6.25"	3/0
104-15556	Webster STILLE Delicate Smooth 12 cm 4.75"	4/0	104-15573-181	Ryder STILLE TC jaws 18 cm 7"	3/0
104-15557	Webster-Kleinert STILLE Smooth Offset bows 12 cm 4.75"	4/0	104-15573-201	Ryder STILLE TC jaws 20 cm 8"	3/0
104-15558	Webster STILLE Smooth TC jaws 12 cm 4.75"	3/0	104-15585-161	STILLE Diamond Jaws 16 cm 6.25"	4/0
104-15559	Webster STILLE Serrated TC jaws 12 cm 4.75"	4/0	104-15585-181	STILLE Diamond Jaws 18.5 cm 7.5"	4/0
104-15560	Halsey STILLE TC jaws 12.5 cm 5"	3/0	104-15585-221	STILLE Diamond Jaws 22 cm 8.75"	4/0
104-15560-2	Halsey STILLE Ergonomic TC jaws 12.5 cm 5"	3/0	104-15585-251	STILLE Diamond Jaws 25 cm 10"	4/0
104-15561-141	Mayo-Hegar STILLE TC jaws 14 cm 5.5"	0	104-15586	STILLE TC jaws 18 cm 7"	4/0
104-15561-161	Mayo-Hegar STILLE TC jaws 16 cm 6.25"	0	104-15590-121	Olsen-Hegar STILLE TC jaws With scissors 12 cm 4.75"	3/0
104-15561-181	Mayo-Hegar STILLE TC jaws 18 cm 7"	0	104-15590-131	Olsen-Hegar STILLE TC jaws With scissors 13.5 cm 5.25"	2/0
104-15561-201	Mayo-Hegar STILLE TC jaws 20 cm 8"	0	104-15590-171	Olsen-Hegar STILLE TC jaws With scissors 17 cm 6.75"	0
104-15561-251	Mayo-Hegar STILLE TC jaws 25 cm 10"	0	104-15590-181	Olsen-Hegar STILLE TC jaws With scissors 18.5 cm 7.25"	0
104-15562-161	Mayo-Hegar STILLE Heavy Pattern TC jaws 16 cm 6.25"	1	104-15621-179	Mathieu STILLE TC jaws 17 cm 6.75"	1
104-15562-181	Mayo-Hegar STILLE Heavy Pattern TC jaws 18 cm 7"	1	104-15621-209	Mathieu STILLE TC jaws 20 cm 8"	1
104-15562-201	Mayo-Hegar STILLE Heavy Pattern TC jaws 20 cm 8"	1	104-15665-151	Frimand-STILLE Duo-Grip Needle holder TC jaws 15 cm 6"	2/0
104-15562-241	Mayo-Hegar STILLE Heavy Pattern TC jaws 24 cm 9.5"	1	104-15777	Kaye STILLE With scissors Serrated TC jaws 12 cm 4.75"	5/0
104-15562-271	Mayo-Hegar STILLE Heavy Pattern TC jaws 27 cm 10.5"	1	104-15785	Converse STILLE TC jaws 12 cm 4.75"	4/0
104-15562-301	Mayo-Hegar STILLE Heavy Pattern TC jaws 30 cm 12"	1			
104-15565-151	Crile-Wood STILLE TC jaws 15 cm 6"	2/0			
104-15565-181	Crile-Wood STILLE TC jaws 18 cm 7"	2/0			
104-15565-201	Crile-Wood STILLE TC jaws 20 cm 8"	2/0			
104-15565-231	Crile-Wood STILLE TC jaws 23 cm 9"	2/0			

Suture size indicates the largest recommended for each needle holder.



104-15542

DERF STILLE NEEDLE HOLDER | Nålförare Derf STILLE

Needle holder with rather short TC jaws, popular in plastic surgery. Intended for 3/0 suture and finer.


ART. NUMBER	LENGTH	TC
104-15542	12 cm 4.75"	•



Art. Number: 104-15542



104-15545



104-15546

WEBSTER STILLE DELICATE NEEDLE HOLDER | Nålförare Webster STILLE gracil modell

Very delicate needle holder for very fine sutures. Smooth jaws or Diamond jaws. Intended for 5/0 suture and finer.

ART. NUMBER	JAWS	LENGTH	DIAMOND
104-15545	Very delicate	12 cm 4.75"	•
104-15546	Very delicate Smooth	12 cm 4.75"	



Art. Number: 104-15545

**KILNER STILLE NEEDLE HOLDER |
Nålförare Kilner STILLE**

Delicate needleholder with bayonet shaped shanks. TC jaws. Intended for 4/0 suture and finer.


ART. NUMBER	LENGTH	TC
104-15549-131	13 cm 5"	•



Art. Number: 104-15549-131

**GILLIES STILLE NEEDLE HOLDER |
Nålförare Gillies STILLE**

Needle holder with scissors. Ergonomic design on shanks and finger rings. One rather large ring handle is set at an angle for fitting the top of the thumb. Without ratchets. With or without TC jaws. The scissor part enables suturing and cutting with one instrument. Intended for 2/0 suture and finer.



104-15550



104-15551


ART. NUMBER	JAWS	LENGTH	TC
104-15550	With scissors	16 cm 6.25"	
104-15551	With scissors	16 cm 6.25"	•



Art. Number: 104-15551

WEBSTER/KLEINERT STILLE NEEDLE HOLDER |
Nålförare Webster/Kleinert STILLE

Delicate needle holder for fine needles. With smooth jaws or TC jaws.
 The Webster-Kleinert (104-15557) has offset bows for a comfortable grip. Intended for 4/0 suture and finer.



104-15556
104-15557

104-15559


ART. NUMBER	TYPE	JAWS	LENGTH	TC
104-15556	Delicate	Smooth	12 cm 4.75"	•
104-15557	Offset bows	Smooth	12 cm 4.75"	•
104-15559	Serrated		12 cm 4.75"	•



Art. Number: 104-15559

WEBSTER STILLE NEEDLE HOLDER |
Nålförare Webster STILLE

Delicate needle holder for fine needles. With smooth TC jaws. Intended for 3/0 suture and finer.



104-15558

ART. NUMBER	JAWS	LENGTH	TC
104-15558	Smooth	12 cm 4.75"	•



Art. Number: 104-15558



104-15560

HALSEY STILLE NEEDLE HOLDER |
Nålförare Halsey STILLE

Needle holder with TC jaws used in many procedures. One model with ergonomically designed finger ring for thumb. For slightly larger needles than the Webster needle holder. Intended for 3/0 suture and finer.

ART. NUMBER	FINGER RING	LENGTH	TC
104-15560		12.5 cm 5"	•
104-15560-2	Ergonomic	12.5 cm 5"	•



104-15561

MAYO-HEGAR STILLE NEEDLE HOLDER |
Nålförare Mayo-Hegar STILLE

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15561-141	14 cm 5.5"	•
104-15561-161	16 cm 6.25"	•
104-15561-181	18 cm 7"	•
104-15561-201	20 cm 8"	•
104-15561-251	25 cm 10"	•



Art. Number: 104-15560



Art. Number: 104-15561-181

MAYO-HEGAR STILLE HEAVY NEEDLE HOLDER | Nålförare Mayo-Hegar STILLE kraftig modell

Heavier version of the widely used Mayo-Hegar needle holder, with TC-jaws. Intended for 1 suture and finer.

ART. NUMBER	PATTERN	LENGTH	TC
104-15562-161	Heavy	16 cm 6.25"	•
104-15562-181	Heavy	18.5 cm 7.25"	•
104-15562-201	Heavy	20.5 cm 8"	•
104-15562-241	Heavy	24 cm 9.5"	•
104-15562-271	Heavy	27 cm 10.75"	•
104-15562-301	Heavy	30.5 cm 12"	•




104-15565


CRILE-WOOD STILLE NEEDLE HOLDER | Nålförare Crile-Wood STILLE

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15565-151	15 cm 6"	•
104-15565-181	18 cm 7"	•
104-15565-201	20 cm 8"	•
104-15565-231	23 cm 9"	•



Art. Number: 104-15562-201



Art. Number: 104-15565-201



104-15566

CRILE-WOOD STILLE DELICATE NEEDLE HOLDER | Nålförare Crile-Wood gracil modell

Delicate version of the popular Crile-Wood needle holder. With TC jaws. Intended for 3/0 suture and finer.

ART. NUMBER	PATTERN	LENGTH	TC
104-15566-151	Delicate	15 cm 6"	•
104-15566-181	Delicate	18 cm 7"	•
104-15566-201	Delicate	20 cm 8"	•
104-15566-231	Delicate	23 cm 9"	•



104-15665-151

STILLE-FRIMAND DUO-GRIP NEEDLE HOLDER | Nålförare STILLE-Frimand Dou-Grip

Designed in cooperation with Carl-Fredrik Frimand, M.D. Intended for precise use in the palmar grip without the need of grip change during suturing. This new and innovative needle holder reduces the suturing time and tissue manipulation as well as offering better ergonomics for the surgeon. Intended for 2/0 suture and finer.


ART. NUMBER	LENGTH	TC
104-15665-151	15 cm 6"	•



Art. Number: 104-15566-151



Art. Number: 104-15665-151



104-15571

NEW ORLEANS STILLE NEEDLE HOLDER | Nålförare New Orleans STILLE

Short tapered TC jaws, one jaw straight serrated and the other cross-serrated. Intended for 3/0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15571	18 cm 7"	•



104-15573

RYDER STILLE NEEDLE HOLDER | Nålförare Ryder STILLE

Needle holder with very narrow, non-tapered, serrated TC jaws. Popular in vascular surgery. Intended for 3/0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15573-161	16 cm 6.25"	•
104-15573-181	18 cm 7"	•
104-15573-201	20 cm 8"	•




Art. Number: 104-15571



Art. Number: 104-15573-181






104-15585

STILLE NEEDLE HOLDER | Nålförare STILLE

Delicate Diamond Jaws. Intended for 4/0 suture and finer.

ART. NUMBER	LENGTH	
104-15585-161	16 cm 6.25"	•
104-15585-181	18.5 cm 7.5"	•
104-15585-221	22 cm 8.75"	•
104-15585-251	25 cm 10"	•




Art. Number: 104-15585-181



104-15586

**STILLE TC NEEDLE HOLDER |
Nålförare STILLE TC**

Delicate needle holder with TC jaws. Intended for 4/0 suture and finer.

ART. NUMBER	LENGTH	
104-15586	18 cm 7"	•



Art. Number: 104-15586


OLSEN-HEGAR STILLE NEEDLE HOLDER | Nålförare Olsen-Hegar STILLE

Popular needle holder with scissors for suturing and cutting with the same instrument. With TC jaws. See table "Needle holder in relation to suture size" for suitable suture size for each article code. 3+3 ratchets, 104-15590-171 has 2+2 ratchets.

ART. NUMBER	JAWS	LENGTH	TC
104-15590-121	Scissors	12 cm 4.75"	•
104-15590-131	Scissors	13.5 cm 5.25"	•
104-15590-171	Scissors	17 cm 6.75"	•
104-15590-181	Scissors	18.5 cm 7.25"	•



Art. Number: 104-15590-131



104-15621

MATHIEU STILLE NEEDLE HOLDER | Nålförare Mathieu STILLE

Needle holder with straight TC jaws. Plier type handle with ratchets that allows the surgeon to open and close the needle holder rapidly. Intended for 1 suture and finer.

ART. NUMBER	LENGTH	TC
104-15621-179	17 cm 6.75"	•
104-15621-209	20 cm 8"	•



Art. Number: 104-15621-179



Art. Number: 104-15777



104-15777

KAYE STILLE NEEDLE HOLDER |
Nålförare Kaye STILLE

Delicate needle holder with short, fine TC jaws and sharp scissors. The scissor part enables suturing and cutting with one instrument. Intended for 5/0 suture and finer.

ART. NUMBER	JAWS	LENGTH	TC
104-15777	Serrated	12 cm 4.75"	•



Art. Number: 104-15785

CONVERSE STILLE NEEDLE HOLDER |
Nålförare Converse STILLE

Short and delicate needle holder for fine work. With TC jaws. Intended for 4/0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15785	12 cm 4.75"	•

I still have STILLE needle holders that date back from the 1970's and they are performing just as well as the new recent ones I purchased.

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon



104-15140

RAPP-STILLE KNOT-HOLDING FORCEPS |
Knuthållartång Rapp-STILLE

A 13 cm knot-holding forceps with smooth jaws. Ideal instrument to use as a tool for holding the suture knot in narrow areas.

ART. NUMBER	LENGTH
104-15140	13 cm 5"




Art. Number: 104-15140

STILLE WIRE TWISTING FORCEPS | Ståltrådsroterare STILLE

A good tool for twisting steel wires. With TC jaws.

ART. NUMBER	LENGTH	TC
104-15832	15 cm 6"	•
104-15835	18.5 cm 7.25"	•




Art. Number: 104-15835

STILLE WIRE HOLDING FORCEPS | Ståltrådshållare STILLE

For holding steel wires. Equipped with double-action joint and TC jaws. The jaws close parallel to secure a firm grip of the wire.

ART. NUMBER	SIZE	LENGTH	TC
104-15854	Fine	18 cm 7"	•
104-15855	Medium	18 cm 7"	•
104-15856	Large	18 cm 7"	•



Art. Number: 104-15854